Wash-down Modular Pump with Wall-mount Controller  
system 77110-40

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## SAFETY PRECAUTIONS

**DANGER:** High voltages exist and are accessible in the Modular Drive. Use extreme caution when servicing internal components.

**WARNING:** Tubing breakage may result in fluid being sprayed from pump. Use appropriate measures to protect operator and equipment.

**WARNING:** To reduce risk of injury, do not pump materials hotter than 150 degrees Fahrenheit, (65.5°C).

**WARNING:** Power must be removed from pump before removing or installing tubing. Fingers or loose clothing could get caught in drive mechanism.

**CAUTION:** To avoid electrical shock, the power cord protective grounding conductor must be connected to ground. Not for operation in wet locations as defined by EN 61010-1.

### Explanation of Symbols

**CAUTION:** Risk of Danger. Consult Operator's manual for nature of hazard and corrective actions.

**CAUTION:** Risk of crushing. Keep fingers away from rotor while pump is in operation. Stop pump before loading or unloading tubing.

**CAUTION:** Hot Surface. Do not touch.

**CAUTION:** Risk of electric shock. Consult Operator's manual for nature of hazard and corrective actions.

### WARNING: PRODUCT USE LIMITATION

This product is not designed for, nor intended for use in patient connected applications; including, but not limited to, medical and dental use, and accordingly has not been submitted for FDA approval.

## Introduction

The Modular Drive controls the speed of MASTERFLEX® B/T® (Batch/Transfer) pumps to provide flow rates from 0.01 to 3,785 mL.

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## Control/Display Functions/ Mounting Dimensions

Press buttons to activate function.

Use up/down (▲, ▼) arrows to correct/change a flashing display.

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**A) DOWN ARROW (DECREMENT)**—Decrease value of a flashing display.

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**B) UP ARROW (INCREMENT)**—Increase value of a flashing display.

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**C) DISPENSE/COPY**—Set dispense volume, copy amount, or dispense time.

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**D) FLOW CONTROL**—Set flow rate for selected tubing size. To change flow rate, press ▲ or ▼ arrows. (If pump is running, its speed will change with new settings.)

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**E) CAL CONTROL**—Refine built-in calibration, using a measured volume.

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
**G) PRIME**—Run pump at full speed to fill or clear lines.

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## CONTROL/DISPLAY FUNCTIONS/MOUNTING DIMENSIONS

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H) **DIRECTION**—To change motor direction.

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I) **MODE SELECT**—INT for internal control; mA for remote current control; V for remote voltage control.

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J) **SIZE**—Select tubing size and flow units, also displays maximum flow rate.

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
**Mounting Dimensions for Model Numbers 77110-40, -47**

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## Tubing Types

Use only MASTERFLEX® PERFECTPOSITION® B/T® tubing with MASTERFLEX pumps to ensure optimum performance. Use of other tubing may void applicable warranties.

**NOTE:** Use MASTERFLEX PERFECTPOSITION B/T tubing. These pumps are designed to use tubing sizes 87 and 91 only. The tubing sizes refer to the last two digits of the MASTERFLEX tubing model number. All MASTERFLEX tubing in these sizes can be used with this pump. Be sure tubing material matches application.

The following chart provides information for the various sizes of tubing.

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Characteristic	Tubing Size	
	87	91
Inside Dia. in (mm)	0.5 (12.7)	0.75 (19.0)
Outside Dia. in (mm)	1.0 (25.4)	1.25 (31.8)
Hose barb size in (mm)	1/2 (12.7)	3/4 (19.0)
Flow Range (with 11–321 rpm drive)	0.6–19.0 liter/min (0.17–5.0 GPM)	1.2–36.0 liter/min (0.3–9.5 GPM)
Nominal Flow Per Revolution	70.46 mL	141 mL
Maximum Vacuum	28.5 in Hg	32 in Hg
Maximum; Outlet Pressure	35 psi (1.4 bar)	30 psi (1.4 bar)

Flow ranges are approximate-calculated under the following conditions: 0 psi at inlet, 0.5 psi at outlet and water circulant at 22°C (72°F).

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## INSTALLING THE PUMP TUBING

**WARNING:** Power must be removed from pump before removing or installing tubing. Fingers or loose clothing could get caught in drive mechanism. Do not operate this pump without cover or interlock door properly closed and latched. Rotating parts can cause serious injury.

1. Cut off power to the pump by disconnecting line cord or, if wired permanently, by removing the fuse. Do not assume that turning off the switch at the motor (or controller) is "safe enough."
2. Unlatch the locking clamp and open the cover.
3. Insert the tube in the appropriate upper tube retaining pocket (see Figure 3). Line up the "PERFECTPOSITION" placement marks printed on the tube with the outside edge of the retainer assembly (see Figure 1).
4. Wrap the tubing around the rotor assembly and insert the tubing in the lower retaining pocket (see Figure 2).
5. Close the cover and re-latch the locking clamp.
6. If the new tube must be cut from a length of approved replacement tubing, a minimum of 32 inches will be required for a new tube.

**WARNING:** Do not operate this pump without cover or interlock door properly closed and latched. Rotating parts can cause serious injury.

7. Restore power to the pump.

**NOTE:** Under some circumstances, tubing may creep into pump. If this problem occurs it can be remedied by installing a hose clamp or fitting immediately upstream of the upper retaining pocket.

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
## REPLACING THE PUMP TUBING REMPLACEMENT DES TUBES

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
**Figure 3.**  
**Figure 3.**  
**Abbildung 3.**  
**Figure 3.**  
**Figure 3.**

A) TUBE RETAINING POCKETS FOR  
(A) 87 AND (B) 91 TUBING  
POCHES DE RETENUE DE TUBE  
POUR TUBE (A) 87 ET (B) 91  
SCHLAUCHHALTERUNG FÜR  
(A) 87 UND (B) 91 SCHLAUCHE  
CAVIDADES DE RETENCIÓN DE  
TUBOS PARA TUBOS (A) 87 Y (B) 91  
TASCHE DEL TUBO PER TUBI DA (A)  
87 E (B) 91

## Replacing Rollers

To replace rollers:

1. Using a retaining ring tool, part no. 07585-06, remove the retaining rings from the ends of the roller axles and slide the rollers off. **Take care to avoid opening the retaining rings too wide.**
2. Check to be sure that the wave washers are installed on the axles against the rotor plate.
3. Slide the new rollers onto the axles, placing the ends with the flush bearing surface inward toward the rotor plate against the wave washers, and the etched ends with the recessed bearing surface outward toward the free ends of the roller axles.
4. Replace the retaining rings. You may have to push the rollers in to compress the wave washers to allow the retaining rings to engage the grooves in the axles.

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
## INTERCHANGING ROLLERS PUMP MOUNTING DIMENSIONS

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
## Controller Setup

1. Connect Motor Cable plug to mating receptacle on the Controller.
2. Connect power cord of Controller to grounded power line outlet.
3. Turn pump on and select TUBING SIZE.

**NOTE:** If CAL LED is lit, that tubing size has been previously field calibrated. If LED is not lit, the drive is operating with the built-in factory calibration. To clear a field calibration, press and hold the CAL switch until the CAL light goes out. This will take about 3 seconds. To recalibrate for better accuracy, see Calibration section.

4. MODE selection (INT, mA, V).
5. Select MOTOR DIRECTION (CW or CCW).
6. PRIME and CALibrate the pump (if required).
7. Press FLOW key and watch display to set the flow rate with UP/DOWN keys.
8. Press STOP/START key to begin pumping.

**NOTE:** Under some circumstances, tubing may creep into pump. If this problem occurs it can be remedied by installing a hose clamp or fitting immediately upstream of and very close to the inlet port.

**NOTE:** Pump will restart automatically after a brownout or powerout condition.

## Calibration

1. Select correct tubing size and flow rate.
2. Press CAL, calibration volume appears.
3. Press STOP/START, the pump will use its stored memory to dispense the specified calibration sample quantity. The pump will stop automatically.
4. Weigh/measure the sample.
5. Use UP/DOWN arrow keys to correct the flashing display.

**NOTE:** If the adjusted calibration is too great, "Err" will appear in the display. If this occurs, press the CAL control and repeat the calibration procedure. The microprocessor will retain one special calibration value per tubing size, even when power is turned off. The next calibration will replace the existing value.

6. Press SIZE to exit the calibration cycle.

## Maximum Flowrate (OTHER Tubing)

1. To set the maximum flowrate for non-standard pumps or tubing sizes, OTHER press CAL, then FLOW. The maximum flowrate will then flash on the display.
2. Use UP/DOWN arrow keys to set desired flowrate.
3. Press SIZE to exit.

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## DISPense/copy

A first press of the DISP key results in the last entered dispense volume being displayed. The "VOL" annunciator will illuminate and flash. The INC/DEC keys are used to change the dispense volume, if desired. The STOP/START key then initiates delivery of the set volume. The amount remaining to be dispensed will be displayed during countdown. The dispense function is exited by pressing any key except Increment, Decrement, DISP, or STOP/START.

A second press of the DISP key causes the COPY annunciator to illuminate and flash. The STOP/START key is then used to set the desired volume without the need to know the volume in specific units. A third press of the DISP key enters the volume dispensed. The COPY annunciator stops flashing. The STOP/START key is then used to initiate delivery of the copied volume. The number of copies dispensed will be displayed after each dispense. The STOP/START key is used to pause the copy dispense during dispensing; copy dispense can then be continued using the STOP/START key.

A fourth press of the DISP key results in the last entered dispense time being displayed. The SEC annunciator will illuminate and flash. The INC/DEC keys are used to change the dispense time, if desired, from 1 to 9999 seconds. The STOP/START key then initiates delivery for the set time interval. The remaining time will be displayed during countdown. Pressing the DISP key a fifth time exits this mode.

## Keypad Lockout Enable/Disable

Press and hold FLOW. After five (5) seconds, display will change to all dashes. Then, while holding FLOW, press PRIME five (5) times.

The MODE annunciator will flash when the keypad is locked.

## REMOTE CONTROL

Selectable input (0–20 mA, 4–20 mA, 0–10V DC )  
±0.5% linearity control

2300V isolation potential

STOP/START; CW/CCW; PRIME via contact closure

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## Remote Control Setup

1. Place the power switch in the off position.

**CAUTION:** Power must be turned off before connecting the external remote control cable to prevent damage to the drive.

2. Connect the cable from the external remote control to the mating receptacle on the bottom panel.
3. Select type of remote control input and output required as follows:
  - a). Press and hold the MODE key while turning the power switch to the on position. After two seconds, release the MODE key. The initial display will show: "inP". After two seconds the display will show either 0-20 or 4-20.

**NOTE:** Press the up (increment) or down (decrement) arrows to select between 4-20 and 0-20 for current loop control.

- b). Press the MODE key again. The initial display will show: "out". After two seconds the display will show either 0-20, 4-20, or 0-10.

**NOTE:** Press the up (increment) or down (decrement) arrows to select between 4-20 and 0-20 for current loop output, or 0-10 for voltage output.

4. Press the MODE key to select mode of operation. The LED's indicate the selected mode. Select either mA or V.

**NOTE:** If only remote STOP/START, PRIME and/or CW/CCW is to be used, the MODE control can be set to any of the three positions.

5. To adjust the voltage or current scaling for other than zero to full scale:
  - a) Press the MODE key and then the FLOW key at the same time. The display will show "LO" and then the flow rate for minimum current/voltage (factory default = 0).
  - b) Use the UP/DOWN arrow keys to change the flow rate for minimum current/voltage.
  - c) Press the FLOW key. The display will show "HI" and then the flow rate for maximum current/voltage. Use the UP/DOWN arrow keys to change the flow rate for maximum current/voltage or press the SIZE key to set it to maximum flow rate (factory default). Press any other key to save and exit.

The same scaling will be used for both input and output. Each tube size has its own scaling.

**NOTE:** The maximum flow rate for a tubing will change after a calibration is performed. To retain control of the entire flow range, the "HI" scale setting must be changed to the new maximum flow rate after a calibration is performed.

6. Remote STOP/START can be configured to be optional ("OFF") or mandatory ("ON"). When "ON" is selected, drive will not run unless remote STOP/START is closed. When "OFF" is selected (factory default), remote STOP/START can be used to start drive, but drive can also be started by keypad or remote inputs when remote STOP/START is open. Internal mode or remote mode (mA or V) each have their own STOP/START configuration, so first select the desired operating mode before changing STOP/START setup.
  - a) Press and hold the MODE key until the display changes to "STOP". The display will alternate with an "ON" or "OFF".
  - b) Use the UP/DOWN arrow keys to select "ON" or "OFF".
  - c) Press any other key to save and exit.

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
### CONFIGURACION DE CONTROL REMOTO

### MESSA A PUNTO DEL TELECOMANDO

## Connector Pin configuration/with wiring scheme

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





**NOTE:** Colors are those of Remote Cable, Cat. number MN-77300-32.

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## TROUBLESHOOTING AND MAINTENANCE

**DANGER:** Modular Drive Systems are not customer serviceable, except for replacement of FUSES, MOTOR BRUSHES and MOTOR.

**DANGER:** High voltages exist and are accessible in the Modular Drive. Refer servicing to qualified personnel.

**WARNING:** Unplug Controller prior to checking/replacing FUSES, MOTOR BRUSHES and MOTOR.

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
### Tools required:

56 C frame mounting, 5/16 Hex key & (2) 1/2" wrenches.

ISO-71 frame mounting, 10 mm wrench & (2) 13 mm wrenches.

To install replacement motor refer to Figure 4 and follow these steps:

1. Apply anti-seize compound to shaft and key. Slide motor forward to engage the male motor coupling with the female gear head coupling. Rotate pump rotor if necessary to align the couplings.
2. Using the 5/16 inch hex key, bolt the motor to the gear head adapter by inserting the four bolts through the gear head adapter into the motor. Torque bolts from 17 to 19 foot-pounds.
3. Using the 1/2 inch wrench, secure the motor mounting feet to the base using the four bolts. Torque bolts from 9 to 11 foot-pounds.

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POMPA

## Motor Brush Replacement

**Tools Required:** Phillips screwdriver

**WARNING:** Power must be removed from motor before performing this procedure

1. Cut off power to the motor by disconnecting line cord. Do not assume that turning off the switch at the motor (or controller) is "safe enough."
2. To access the motor brushes, remove six screws securing the access plate and remove plate (see Figure 5).
3. Loosen screw terminal at top of brush housing and disconnect brush wire.
4. Press down on brush retainer to disengage tabs then rotate brush retainer slightly toward front of motor and remove brush retainer.
5. Slide brush assembly out of housing.
6. Install new brush assembly with brush wire toward rear of motor and spring assembly on top.
7. Insert brush retainer against brush springs and push down, then rotate retainer slightly toward rear of motor to engage tabs of retainer under rear edge of brush housing.
8. Attach brush wire to screw terminal at top of brush housing. Be sure wire is clear of access opening.
9. Attach cover plate with six screws.

**NOTE:** Always replace both brushes at the same time.

**NOTE:** Brushes should be checked every 6 months or 2000 operating hours or if erratic operation occurs.

Erratic operation may occur immediately after brush replacement. Allow motor to run up to an hour to allow brushes to seat.

**NOTE:** Replace both brushes, if either brush is less than 7.6 mm (0.300 in) long from base to point.

Contact your dealer if you have service needs.

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## TROUBLESHOOTING AND MAINTENANCE (cont.)

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## Troubleshooting (Cont.)

If an **error message** is displayed, refer to the following list for possible corrective action you can take. If these do not correct the problem, contact your dealer.

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## REPLACEMENT PARTS AND ACCESSORIES

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### SPECIFICATIONS

<b>Output:</b>	
Speed:	11 to 321 r/min
Torque output, Maximum:	1440 oz-in (104 kg•cm)
<b>Speed regulation:</b>	
Line	±0.25% F.S.
Load	±0.25% F.S.
Drift	±0.25% F.S.
Display:	Four-digit, seven segment LED
<b>Remote outputs:</b>	
	Voltage speed output (0–10V DC)
	Current speed output (0–20 mA or 4–20 mA)
	Tach output (TTL, 128 to 7680 Hz)
	Motor running output (N.O. & N.C. contact closure)
<b>Input:</b>	
<b>Supply voltage limits:</b>	
77110-40	90 to 130 Vrms @ 50/60 Hz
77110-47	200 to 260 Vrms @ 50 Hz
<b>Current, max.:</b>	
77110-40	6.5A
77110-47	3.3A
<b>Remote Inputs:</b>	
	START/STOP, CW/CCW, PRIME (contact closure)
	Voltage input (0–10V DC)
	Current input (0–20 mA or 4–20 mA)
<b>Construction:</b>	
<b>Dimensions (L × W × H):</b>	
940-2611, 940-2613	9 in × 11 in × 4 1/2 in (229 mm × 279 mm × 114 mm)
940-2610, 940-2615	26 in × 13 in × 10 5/8 in (660 mm × 330 mm × 270 mm)
<b>Weight:</b>	
940-2611, 940-2613	9.4 lbs (4.3 kg)
940-2610, 940-2615	89.0 lbs (40.37 kg)
<b>Enclosure Rating:</b>	
940-2611, 940-2613	IP 56 (NEMA 4) Per IEC 60529
<b>Environment:</b>	
Temperature, Operating:	32° to 104°F (0° to 40°C)
Temperature, Storage:	–49° to 149°F (–45° to 65°C)
Humidity (non-condensing):	10% to 90%
Altitude:	Less than 2000 m
Pollution Degree:	Pollution Degree 3 (Sheltered locations)
<b>Chemical Resistance:</b>	
	Exposed material is painted aluminum, plastic and vinyl
<b>Compliance:</b>	
	115V: UL508C
	CSA C22.2 No. 14
	230V (for CE Mark):
	EU Low Voltage Directive (EN61010-1) and
	EU EMC Directive (EN61326)

### Cleaning

Keep the drive enclosure clean with mild detergents.  
Never immerse nor use excessive fluid.

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NOTES

## WARRANTY

Use only MASTERFLEX PERFECTPOSITION precision tubing with MASTERFLEX pumps to ensure optimum performance. Use of other tubing may void applicable warranties.

The Manufacturer warrants this product to be free from significant deviations from published specifications. If repair or adjustment is necessary within the warranty period, the problem will be corrected at no charge if it is not due to misuse or abuse on your part, as determined by the Manufacturer. Repair costs outside the warranty period, or those resulting from product misuse or abuse, may be invoiced to you.

The warranty period for this product is noted on the Warranty Card.

## PRODUCT RETURN

To limit charges and delays, contact the seller or Manufacturer for authorization and shipping instructions before returning the product, either within or outside of the warranty period. When returning the product, please state the reason for the return. For your protection, pack the product carefully and insure it against possible damage or loss. Any damages resulting from improper packaging are your responsibility.

## TECHNICAL ASSISTANCE

If you have any questions about the use of this product, contact the Manufacturer or authorized seller.

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Vernon Hills, Illinois U.S.A. 60061-1844